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1. Specifications and Scope of Supply

| | |
|------------------|-------------------------------------|
| Item no.: | PM-2090 |
| Dimensions: | 160 x 60 cm |
| ILN: | 4250355320902 |
| Item no.: | PM-2095 L |
| Dimensions: | 175 x 60 cm |
| ILN: | 4250355320957 |
| Item no.: | PM-2097 XL |
| Dimensions: | 185 x 70 cm |
| ILN: | 4250355320971 |
| Material: | PES/Carbon, latex free |
| Color: | blue/red |
| Storage: | -5 – 40° C protect from sunlight |

Scope of supply:

- Caremaster sliding mat
- Manual

Made in Germany

max.
200 kg

2. Use / Placement

The Caremaster sliding mat is used for transfers for patients in seated and supine positions and for turning patients in bed. The Caremaster with the three extra-strong webbed loops can also be used for very heavy persons. The slide mat has a wide range of uses. Please find some examples of use below:

Application 1:

The loops should be placed in such a way that you can easily pull the loops in the required direction. Place the patient in the Sims' position. Fold the half side of the Caremaster each lengthwise and widthwise and place the folded part as close as possible to the back of the patient. Position the patient on the reverse side then by placing the patient from the supine position to the opposite Sims' position and unfold the folded Caremaster then under the patient. Transfer the patient by pulling the belt loops from the bed's head end or from the lateral position by the caregivers. The patient can support the transfer himself by slightly pushing off with the feet (=> help people to help themselves)

Application 2:

The loops should be placed in such a way that you can easily pull the loops in the required direction. Place the patient in the Sims position. Then fold the half side of the Caremaster lengthwise and place the folded part as close as possible to the back of the patient. Position the patient on the reverse side then by placing the patient from the supine position to the opposite Sims' position and unfold the folded Caremaster then under the patient. Position the patient again in a supine position. His or her back is placed on the folded Caremaster now. Transfer the patient by pulling the belt loops from the bed's head end or from the lateral position by the caregivers.

Application 3:

The situation can be managed even more easily with the Alpha® sliding board PM-5050 particularly suitable for bed – wheelchair / wheelchair – car transfer. Fold the Caremaster lengthwise around the slide board. Use the slide board to assist with the transfer as much as possible and bring the Caremaster under the seated patient as far as possible with light pressure. Pull the slide board out from the Caremaster. The patient is now seated on the folded Caremaster. Push the slide board as far as possible under the folded Caremaster and the seated patient. Transfer the patient by pulling the belt loops.

This product is suitable for re-use.

3. Maintenance / Care

The product does not require any special maintenance. It can be washed in a washing machine at 90° C and is suitable for low dryer temperatures (45 – 50° C). When using non-standard washing agents, bleaching agents, softener or other chemical washing additives, the functional properties of the product can strongly degrade. The robust material can be disinfected with the usual wipe disinfectants. Check the compatibility before use where required on the basis of the disinfectant's instructions. We recommend the use of alcohol- and aldehyde-free agents, e.g. *Trionic®*. There is also some care information on the label.

4. Warranty

- Petermann gives a warranty that the product is devoid of material defects and manufacturing defects.
- The warranty is valid for a period of **3 years**. It does not cover any defects which are due to improper use.
- Petermann cannot be held liable for any accidental or indirect injury or damage.

5. Visual and Function Inspection / Safety Information / Re-use

Repositioning products are also subject to natural aging. Check the equipment regularly and carefully to make sure it is working properly and especially before re-use. Make sure that the product is in good condition and that the material and the seams are not damaged. Check the function of the product by applying pressure and loading to the handles and make sure that the material resists the load. The smooth surface of the low-friction material makes it easy to slide and move. The best way to test these properties is by comparing the product with a new product.

For re-using the product it must be washed and disinfected before, see 3.

WARNING:

- ***Faulty products must no longer be used.***
- ***If the patient's mobility is restricted, it is essential to supervise the motions of transfer by providing trusted help.***
- ***The function of the product may be adversely affected if it is washed with non-standard detergents or bleaches, fabric softeners or other chemical detergent additives.***
- ***The smooth surface considerably minimises resistance to friction during the change of position. Carefully check the transfer surfaces, making sure they are stable and brakes are applied (e. g. wheelchair, bed).***

6. Product Labeling

The Petermann Alpha® open sliding mat Caremaster is a medical product, class I and complies with the MDD directive 93/42/EEC.

Petermann GmbH is certified in accordance with the directive DIN EN ISO 13485 by the German TUV.

7. Disposal Information

For the Caremaster there are no special disposal rules; it can be disposed of in accordance with national and local regulations as necessary.

If you have any questions about the use of our products, please call our service number:

+49 (0) 9868 9339-0

Petermann GmbH • Ziegelhaus 12 • D-91601 Dombühl
Postal address: Petermann GmbH • Waldweg 3 • D-91601 Dombühl
Phone: +49 (0) 9868 9339-0 • Fax: +49 (0) 9868 9339-20
E-mail: petermann@pm-med.de • www.pm-med.de

"Information subject to technical modifications to our products and deviations in colour or printing errors"

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Petermann GmbH est certifié par le TUV conforme à la directive DIN EN ISO 13485.
Adresse postale: Petermann GmbH • Waldweg 3 • D-91601 Dombühl